





A device that patients and practices will love

SureSmile® VPro™ is a high-frequency vibration device that patients can use at home to help aligners fit properly. It may reduce discomfort and may reduce treatment time when used five minutes a day with clear aligners.^{1,2}


High-Frequency Vibration Clinical Data

HFV Clinical Overview

When used correctly, SureSmile® VPro™ may reduce discomfort[†]



SureSmile VPro may improve aligner tracking^{†*}



[†]Over not using the device.

Reference: 1. Alansari S, Atique MI, Gomez JP, et al. The effects of brief daily vibration on clear aligner orthodontic treatment. J World Fed Ortho. 2018;7(4):134-140. doi:10.1016/j.ejwf.2018.10.002

2. Shipley T, Farouk K, El-Bialy T. Effect of high-frequency vibration on orthodontic tooth movement and bone density. J Orthod Sci. 2019;8:15. doi:10.4103/jos.JOS_17_19

Intended for use by the orthodontic patient during treatment with aligners to facilitate minor anterior tooth movement.

Happy Patients, Healthy Smiles

SureSmile[®] VPro[™]

Case Examples

CASE No. 1

Diagnosis

- Deep bite
- Deviated midline
- Mandibular crowding and rotations
- Lingual rolling of mandibular molars

Modality

38 clear aligners &
SureSmile VPro

Treatment Time

Expected: 19 months*
Actual: 6 months



Case courtesy of Dr Bill Dischinger.
Oregon, USA

CASE No. 2

Diagnosis

- Class I
- Crowding

Modality

32 clear aligners &
SureSmile VPro

Treatment Time

Expected: 12 to 14 months†
Actual: 5 months



Case courtesy of Dr David R. Boschken.
Oregon, USA

Doctors may have been compensated for providing the above case examples.

*Expected treatment time calculated based on 2-week aligner interval changes.

Expected treatment time provided by treating clinician.

Patients may have been treated with an earlier VPro model.

Recapture During Active Treatment and Retention

Overcoming Non-Compliance

It's difficult for any patient to be 100% compliant with their aligner wear. But the SureSmile® VPro™ Device, when used at home for 5 minutes a day, can help overcome tracking issues.¹

Download the SureSmile VPro App

Patients can pair the SureSmile VPro app with their device in order to:

- Track daily SureSmile VPro usage
- Monitor SureSmile VPro battery life
- Receive reminders to use device daily and change aligners as prescribed
- Share results and progress photos



CASE EXAMPLE

Recapture from relapse during retention

Problem: Relapse due to non-compliant retainer wear

- Costly retreatment may have lost the orthodontic referral of the patient's family member

Solution: Retainers and SureSmile VPro

- Patient reported relief from discomfort
- Seating was improved same day
- Patient returned 4 months later with well seated retainers; teeth were back to original finish position



Retainer not seating



Same day seating with device



4 months with device

*Cases courtesy of Shipley Education Services, LLC.
Doctor may have been compensated for providing the above case example.
Patient was treated with an earlier VPro model.

The excellent clear aligner treatment experience

Patient-inspired solutions that will make them smile

SureSmile® helps you keep patients happy from their first tray to their final smile. SureSmile® Aligners provide predictable results that give patients confidence in the treatment process. And when patients use SureSmile® VPro™ just five minutes a day with their aligners, they may reduce discomfort and may reduce treatment time.¹⁻³



SureSmile® VPro™:

May provide added comfort^{1,2}

Patients reported less discomfort with the SureSmile VPro device

May reduce aligner treatment time^{1,3}

Using SureSmile VPro may speed up treatment time by up to 50%

May reduce office visits¹

Using SureSmile VPro can help patients complete their treatment with fewer trips to the dentist

Brochure not approved for use in Japan

References: 1. Alansari S, Atique MI, Gomez JP, et al. The effects of brief daily vibration on clear aligner orthodontic treatment. J World Fed Ortho. 2018;7(4):134-140. doi:10.1016/j.ejwf.2018.10.002 2. Shipley T. Effects of high frequency acceleration device on aligner treatment—a pilot study. Dent J (Basel). 2018;6(3):32. doi:10.3390/dj6030032 3. Nicozisis J, Brigham G, Sparaga J, Shipley T. Effects of 120Hz high-frequency acceleration device on orthodontic discomfort. Preprints. Preprint posted online October 24, 2018. doi:10.20944/preprints201810.0557.v1

Contact your Dentsply Sirona rep by calling **1.888.672.6387**
or visit **www.dentsplysirona.com/suresmile** to learn how you
can offer the confidence of an improved smile to your patients.



THE DENTAL
SOLUTIONS
COMPANY™

